


Urinary Sheath Sizing Guide



To ensure the most accurate measurements, this guide must be printed on A4 paper. Any scaling or 'fit to page' options should be deselected. All printers are different, and so this guide should only be used as a sizing estimate.

- Cut out the measuring guide following the lines carefully
- Pass the guide around a flacid penis near the base
- Note down which size the arrow indicates
- Try the sheath that fits the indicated size

The urinary sheath should be a comfortable fit, not too tight as to constrict the penis but not too loose that it could come off.

We recommend trying the marked size first and the next size up before making a considerable purchase to ensure a perfect fit.

The sheaths feature an adhesive inside to help keep it in place. Removal of the sheath should be relatively easy but can be aided with soap and warm water if needed.

Size 1
24mm

Size 2
28mm

Size 3
31mm

Size 4
35mm

Size 5
40mm